

Vital Signs
a GE Healthcare Company

ACTAR*



ACTAR D-fib* Adult Training Manikins

ACTAR D-fib provides one cost-effective manikin for both cardio-pulmonary resuscitation (CPR) and automated external defibrillation (AED) training. Equipped with anatomical landmarks, including an Adam's Apple, carotid groove, nipples and lower rib for hand alignment, the student will practice compressions, pulse check and AED pad placement. The jaw thrust, head tilt, and navel will teach the student proper technique during choking scenarios.

The ACTAR D-fib is available in five, ten and twenty five packs (pictured right) with disposable lungs and carry cases.



ACTAR AED Training System

The AED Trainer is a training kit that contains:

- AED Trainer – utilizes three buttons (green, yellow and red) to simulate on/off, analyze and shock.
- Training Electrodes – adhesive pads that are secured to the chest of the ACTAR manikin.
- Cables – attach directly to the AED trainer and connect to the training electrodes.
- Audio CD – takes the student through various emergency scenarios.



All components come packed inside the padded case and are available in sets of five and ten.



Armstrong Medical is the exclusive dealer for ACTAR products within the United States. For sales outside the United States, contact your local representative.

ACTAR 911*

The basic but effective CPR training manikin comes standard as an adult size, but easily converts to child size with the Cadet Pac (sold separately).

ACTAR Patrol*

Includes five adult manikins, disposable lungs, airways, instructions and video, all packed in an easy to carry nylon bag.

ACTAR Squadron* (pictured)

Includes ten adult manikins, disposable lungs, airways, instructions and video, all packed in an easy to carry nylon bag.

Cadet Pac*

Converts ACTAR 911 manikins into ACTAR 911 child simulators. Includes 5 child-size compression pistons, instruction card and nylon tube. The carrying tube clips onto Squadron or Patrol bags.

Other available accessories include replacement disposable lungs, chest plates, adult pistons and replacement heads. (pictured)



ACTAR Infantry*



Infant CPR manikins available in sets of five or ten, packaged in a convenient nylon carry bag.

Also included are disposable lungs, instructions, and training video. Replacement accessories for the Infantry are also available.

ABOUT GE HEALTHCARE:

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our “healthymagination” vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world. Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

Vital Signs, Inc.
a GE Healthcare Company
20 Campus Road
Totowa, NJ 07512
USA



imagination at work

© 2011 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram, Vital Signs are trademarks of General Electric Company.

* ACTAR, ACTAR Infantry, ACTAR D-fib, ACTAR 911 Squadron, ACTAR 911 Patrol, ACTAR 911 Cadet Pac, ACTAR 911 Trooper, and Actar 911 are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.

L2075_ACTAR_Global_letter_v1_October 2011
DOC1041988
ACCC-0336-10.11-EN-US