

Vital Signs
a GE Healthcare Company

Disposa-View*
Prone Position Head Cushion with Mirror



View, Access and Support intubated patients in the prone position with the single-use Disposa-View head positioner.



- Gain clear view and access to face, eyes, and endotracheal tube through slotted mirror design
- Control infections and eliminate reprocessing with single-use
- Comfort and support patient's head and neck in the contoured face cradle

Product Number	8000HDP
Color	Natural
Dimensions	Height: 5.84 in (14.83 cm) Width: 9.50 in (24.13 cm) Length: 12.05 in (30.61 cm)
Materials	Cushion Material: Polyurethane Foam Mirror material: Acrylic
Tensile Strength	Greater than 15 lbs /sq. in.

Vital Signs, Inc
a GE Healthcare Company
20 Campus Road
Totowa, New Jersey 07512
USA

© 2011 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

*Disposa-View is a trademark of General Electric Company
GE, GE Monogram, and Vital Signs, are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.



imagination at work

L2079_Disposa-View_Global_v1_October 2011
DOC1041997
ACCC-0335-10.11-EN-US