

Vital Signs  
a GE Healthcare Company

## Disposable Face Masks A-Series Masks



## A-Series Masks: Suit your individual preference with the several styles of face masks that our traditional line offers.

Part Number	Description	Case Quantity	Size	Weights (g)
5210 5211	Clear, rigid dome, air-filled cushion face mask with no tail valve	20 100	1	12.2 to 14.8
5220 5221EU	Clear, rigid dome, air-filled cushion face mask with no tail valve	20 100	2	17.2 to 20.1
5230 5231EU	Clear, rigid dome, air-filled cushion face mask with tail valve	20 100	3	24.0 to 27.3
5230X 5231X	Clear, rigid dome, air-filled cushion face mask with no tail valve	20 100	3	21.8 to 25.3
5240 5241	Clear, rigid dome, air-filled cushion face mask with tail valve	20 100	4	32.1 to 39.8
5240X 5241X	Clear, rigid dome, air-filled cushion face mask with no tail valve	20 100	4	30.1 to 37.8
5245 5246	Elongated shape, clear, rigid dome, air-filled cushion face mask with tail valve	20 100	5	49.2 to 57.1
5250EU 5251	Round shape, clear, rigid dome, air-filled cushion face mask with tail valve	20 100	5	46.5 to 55.1

### Materials

Cushion	PVC Plastisol
Face Plate	Polycarbonate or equivalent

### Physical Characteristics

Color	Yellow tint
Tail Valve	5230, 5231, 5240, 5241, 5245, 5246, 5250 and 5251
No Valve	5210, 5211 5220, 5221, 5230X, 5231X, 5240X and 5241X

### Connectors

Face Plate Connector	ISO 15 mm Male (ISO STD 5356-2 2004 (E), size 5210 ISO 22 mm Female (ISO STD 5356-2 2004 (E) size: 5220, 5230, 5230X, 5240, 5240X, 5245 and 5250
Prong Rings	Not included on size 5210, 5220 Attached to 22 mm Female patient connection One on each: 5230, 5230X, 5240, 5240X, 5245 and 5250
Tail Valve	Luer Connector (ISO STD 594)

Vital Signs, Inc.  
a GE Healthcare Company  
20 Campus Road  
Totowa, New Jersey 07512  
USA

© 2012 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram, Vital Signs are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.



imagination at work