

Vital Signs
a GE Healthcare Company

Disposable Face Masks C-Series Masks



C-Series Masks: See the difference in our newest line of face masks featuring color coded hook rings, brilliant clarity, and semi-soft rounded faceplates

Part Number	Size	Prong Color	Connector	Weights (g)	Case Quantity
4810VS	1	Red	15 mm Male	8	20
4820VS	2	Green	15mm Male	11	20
4830VS	3	Purple	22 mm Female	16	20
4840VS	4	Yellow	22 mm Female	20	20
4851VS	5	Blue	22 mm Female	29	100

Materials

Cushion	Polyvinyl Chloride
Face Plate	Polycarbonate
Prong Ring	Polypropylene
Valve	Silicone
Top Valve	Luer Connector (ISO STD 594)

Physical Characteristics

Mask Color	Blow molded clear cushion and face plate
Top Valve	All sizes have a top valve with blue silicone stem

Vital Signs, Inc.
a GE Healthcare Company
20 Campus Road
Totowa, New Jersey 07512
USA

© 2012 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram, Vital Signs are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.

L2077_C Series Face mask_Spec Sheet_Global_May 2012
DOC1150270



imagination at work