

Vital Signs
a GE Healthcare Company

Disposable Face Masks F-Series Flex Masks



F-Series Masks: Feel the quality of our flexible face masks with their soft faceplate, teardrop shape and superior grip.

Part Number	Description	Case Quantity	Size	Weights (g)
6800EU	Clear, flexible dome, air-filled cushion face mask with no top valve	20	0	10.7 to 14.0
6810EU	Clear, flexible dome, air-filled cushion face mask with top valve	20	1	11.4 to 18.4
6820 6821	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	2	12.9 to 22.3
6820BG 6821BG	Bubble gum scent, clear, flexible dome, air-filled cushion face mask with top valve	20 100	2	12.9 to 22.3
6830 6831	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	3	28.0 to 37.0
6830BG 6831BG	Bubble gum scent, clear, flexible dome, air-filled cushion face mask with top valve	20 100	3	28.0 to 37.0
6840 6841	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	4	39.0 to 56.8
6840BG 6841BG	Bubble gum scent, clear, flexible dome, air-filled cushion face mask with top valve	20 100	4	39.0 to 56.8
6850 6851	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	5	45.0 to 67.7
6860 6861	Clear, flexible dome, air-filled cushion face mask with top valve	20 100	6	62.0 to 78.5

Materials

Cushion	PVC Plastisol
Face Plate	Polycarbonate
Scent Additive	Bubble gum scent: 6820BG, 6830BG and 6840BG

Physical Characteristics

Color	Yellow tint
Top Valve	Size: 1, 2, 3, 4, 5 and 6
No Valve	Size: 0

Connectors

Face Plate Connector	ISO 15mm Male (ISO STD 5356-2 2004 (E)) One on each size: 0, 1 and 2 ISO 22mm Female (ISO STD 5356-2 2004 (E)) One on each size: 3,4,5 and 6
Prong Rings	Attached to 22mm Female patient connector One on each: size: 3,4,5 and 6
Top Valve	Luer Connector, (ISO STD 594)

Environmental Conditions

Operational	Temperature: 0 ° F (-18 ° C) to 122 ° F (50 ° C); Humidity: 40 to 95 percent
Storage	Temperature: -40 ° F (-40 ° C) to 140 ° F (60 ° C); Humidity: 40 to 95 percent

Vital Signs, Inc.
a GE Healthcare Company
20 Campus Road
Totowa, New Jersey 07512
USA



imagination at work

© 2012 General Electric Company - All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram, Vital Signs are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.

L2078_F series Face mask_Spec Sheet_Global_v1_May 2012
DOC1150273