



B/V Filter Product Codes	303	303HEPA	5096	5096HEPA
BFE %	99.997	99.9999	99.997	99.9999
VFE %	99.94	99.999	99.94	99.999
Deadspace (ml)	30.0	30.0	27.0	27.0
Weight (g)	16.0	16.0	17.3	18.5
Resistance (cm H <sub>2</sub> O @ 60 LPM)	1.8	2.2	1.8	2.2
Connectors	22 mm F x 22 mm M/15 mm F	22 mm F x 22 mm M/15 mm F	22 mm F/ 15 mm M x 22 mm M/15 mm F	22 mm F/ 15 mm M x 22 mm M/15 mm F
Sampling Port	N/A	N/A	•	•

HCH/F Product Codes	5708	5708HEPA
BFE %	99.997	99.9999
VFE %	99.94	99.999
Deadspace (ml)	51.2	51.2
Weight (g)	20.7	21.8
Resistance (cm H <sub>2</sub> O @ 60 LPM)	2.0	2.4
Moisture Output †	30.5	30.5
Connectors	22 mm F/ 15 mm M x 22 mm M/15 mm F	22 mm F/ 15 mm M x 22 mm M/15 mm F
Sampling Port	•	•

HCH Product Codes	5701	5702	5704	5706
Description	Adult/Pediatric	Adult	Neonatal HCH	Tracheostomy
Deadspace (ml)	30.0	46.0	<5.0	14.0
Weight (g)	13.0	18.0	2.0	6.0
Resistance (cm H <sub>2</sub> O @ LPM)	<1.5 cm H <sub>2</sub> O @ 60 LPM	1.9 cm H <sub>2</sub> O @ 60 LPM	0.2 cm H <sub>2</sub> O @ 2.5 LPM	0.5 cm H <sub>2</sub> O @ 60 LPM
Moisture Output †	29.2	33.2	30.0 † †	18.0
Connectors	22 mm M/15 mm F x 15 mm M	22 mm M/15 mm F x 15 mm M	15 F x 15 M	15 F

† at 500 ml V<sub>T</sub>, 15 BPM in mgH<sub>2</sub>O/L  
† † at 20 ml V<sub>T</sub>

Vital Signs, Inc.  
a GE Healthcare Company  
20 Campus Road  
Totowa, New Jersey 07512 USA



imagination at work

© 2012 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram, Vital Signs are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.

L2087\_VSD Filter Spec Sheet\_Global\_v1\_January 2012

ACCC-0341-01.12-EN-US  
DOC1086761