

Villard

DESIGNER AND MANUFACTURER OF HOSPITAL FURNITURE

eolix 600 x 400



QUALITY REQUIREMENT

www.villard-medical.com


EMERGENCY CART

eolis 600 x 400

Product advantages:

- Manoeuvrability
- Smooth surface (no sharp edges)
- Ergonomy in accordance with standard NFX 35-104 and NFX 35-105
- Compact

Strong and manoeuvrable for quick first emergency situations. Full of external and internal accessories such as cardiac board, central locking system with seals, defibrillator holder (rotating 360°), infusion rod, gaz bottler holder, separation kit drawers...



CRASH CART ref. **1006.01 E01** Exactly equipped as per picture:

- Eolis 600 x 400 (ABS and polymer)
- 2 adjustable side rails (1 left and 1 right)
- Wheels 125 mm + 2 brakes
- Central seal lock (easy cutting)
- 1 small drawer (ht 63 mm)
- 3 medium drawers (ht 138 mm)
- 1 big drawer (ht 216 mm)
- 1 defibrillator tray (360° rotating) + straps
- 1 IV rod
- 1 Cardiac board
- Gaz bottle holder
- Sticker EMERGENCY
- 77 W x 51 D x 101 H

* suggestion of drawer composition

EXTRA ACCESSORIES

Red seals (qty 100) **1100.47**



Holder for suction unit **1103.88**



5 sockets + switch **1101.50**



Other plugs available

Colour choice

Plum

Green

Saffron

Coral

Lilac

Blue

Red

ANAESTHESIA CART

eolis 600 x 400

Product advantages:

- Easy to clean (infection control)
- Storage access

Eolis 600 x 400 cart perfectly designed for anaesthesia. All the accessories have been created to optimize the anaesthesia procedure. The finishing and the material (ABS and Polypropylene) used are especially designed for optimal hygiene and infection control.



ANAESTHESIA CART ref. **1006.01 A01** Exactly equipped as per picture:

- Eolis 600 x 400 (ABS and polymer)
- 3 adjustable side rails (2 left, 1 right)
- Central key lock
- 2 small drawers (ht 63 mm)
- 1 medium drawer (ht 138 mm)
- 2 big drawers (ht 216 mm)
- 1 bridge + medical rail + top shelf
- Sticker ANAESTHESIA
- 77 W x 51 D x 101 H
- Wheels 125 mm + 2 brakes
- Height of the trolley with bridge: 168 cm



ANAESTHESIA CART ref. **1006.01 A02** Exactly equipped as per picture:

- Eolis 600 x 400 (ABS and polymer)
- 3 adjustable side rails (2 left, 1 right)
- Central key lock
- 2 small drawers (ht 63 mm) + separation 12 cases.
- 1 medium drawer (ht 138 mm)
- 2 big drawers (ht 216 mm)
- 1 bridge + 1 medical rail
- Line of 4 boxes + basket
- 1 bin
- 1 sharp container holder
- Catheter holder
- Sticker ANAESTHESIA
- 77 W x 51 D x 101 H
- Wheels 125 mm + 2 brakes
- Height of the cart with bridge: 168 cm

DRAWER KITS

1 level



Kit 2 cases	1600.70	Kit 8 cases	1600.74
Kit 3 cases	1600.71	Kit 12 cases	1600.75
Kit 4 cases	1600.72	Kit 18 cases	1600.76
Kit 6 cases	1600.73	Kit 24 cases	1600.77

2 levels



Kit 2 cases	1600.80	Kit 8 cases	1600.84
Kit 3 cases	1600.81	Kit 12 cases	1600.85
Kit 4 cases	1600.82	Kit 18 cases	1600.86
Kit 6 cases	1600.83	Kit 24 cases	1600.87

3 levels



Kit 2 cases	1600.90	Kit 8 cases	1600.94
Kit 3 cases	1600.91	Kit 12 cases	1600.95
Kit 4 cases	1600.92	Kit 18 cases	1600.96
Kit 6 cases	1600.93	Kit 24 cases	1600.97

SEE MORE ACCESSORIES PAGES 12 TO 15

TREATMENT CART

eolis 600 x 400

Product advantages:

- The accessories can be moved without using tools.
- No metal sliding rails for best hygiene

The Eolis 600 x 400 treatment carts are available in ABS/POLYPROPYLENE. The ISO standard size 600x400 offers you a great volume of storage. Resistant to disinfectants currently used in hospital environment.

TREATMENT CART FULLY EQUIPPED



Ref. **1006.01 T01** Exactly equipped as per picture:

- Cart (ABS + Polymer) with 5 drawers + 2 adjustable side rails (1 left, 1 right)
- Partitions in two drawers (small and medium)
- 1 bin + lateral basket + 3 bottles holder
- 1 sharp container holder
- Moulded worktop + central locking system with key
- Wheels 125 mm + 2 brakes, stainless steel balls bearing
- Sticker TREATMENT
- 77 W x 51 D x 101 H cm

PROCEDURE CART GENERAL CART ISOLATION CART



Ref. **1006.01 G01** Exactly equipped as per picture:

- Cart ABS polymer with 5 drawers + 2 adjustable side rails (1 left, 1 right)
- Adjustable push handle
- Rim for optimal confort
- A large range of accessories is also available to adapt the configuration (procedure, isolation...)
- Wheels 125 mm + 2 brakes, stainless steel balls bearing
- Sticker TREATMENT
- 77 W x 51 D x 101 H cm

OTHER LOCKS



Key lock

Ref. **1600.01**



Seal lock

Ref. **1600.02**



Code lock

Ref. **1600.03**



Badge lock

Ref. **1600.04**

SEE MORE ACCESSORIES PAGES 12 TO 15



Plum



Green



Saffron



Coral



Lilac



Blue



Red

CREATE YOUR CART EOLIS 600 x 400 WITH DRAWERS

eolis 600 x 400


1 THE BASE



10 levels

Ref. **1006.01** - Capacity max. 10 levels

- Lateral and rear colour ABS panels.
- ABS upper beige tray with rims.
- Entirely modular conception offering big flexibility of layout.
- Integrated height-adjustable push handle.
- 2 height-adjustable normalized rails (1 on each side).
- 4 swivel Ø 125 mm polypropylene with threadguards wheels, 2 with brakes, stainless steel ball bearing.
- All around ABS bumper.
- Weight of an empty cart : 27.5 kg



2 LOCKING



Key lock - ref. **1600.01**



Seal lock - ref. **1600.02**



Code lock - ref. **1600.03**



Badge lock - ref. **1600.04**

3 DRAWERS (max. 10 levels)



1 level - ref. **1600.30**
600 x 400 x 63 mm



2 levels - ref. **1600.31**
600 x 400 x 138 mm



3 levels - ref. **1600.32**
600 x 400 x 216 mm

4 SEPARATION KITS



SEPARATION KIT for small drawer 1 level

Kit 2 cases	1600.70	Kit 8 cases	1600.74
Kit 3 cases	1600.71	Kit 12 cases	1600.75
Kit 4 cases	1600.72	Kit 18 cases	1600.76
Kit 6 cases	1600.73	Kit 24 cases	1600.77



SEPARATION KIT for medium drawer 2 levels

Kit 2 cases	1600.80	Kit 8 cases	1600.84
Kit 3 cases	1600.81	Kit 12 cases	1600.85
Kit 4 cases	1600.82	Kit 18 cases	1600.86
Kit 6 cases	1600.83	Kit 24 cases	1600.87



SEPARATION KIT for big drawer 3 levels

Kit 2 cases	1600.90	Kit 8 cases	1600.94
Kit 3 cases	1600.91	Kit 12 cases	1600.95
Kit 4 cases	1600.92	Kit 18 cases	1600.96
Kit 6 cases	1600.93	Kit 24 cases	1600.97

5 COLOURS



Colour choice



Plum



Green



Saffron



Coral



Lilac



Blue



Red

6 ADD EXTRA ACCESSORIES (see pages 12 to 15)

MOBILE WORKSTATION WITH DRAWERS

eolis 600 x 400



Epoxy steel welded monobloc body frame. Closed on 5 sides - Opening side: 600 mm. High-density compact resin work-top (10 mm thick) **with antibacterial treatment**, 99,9% efficiency tested on 7 bacteria according to the norm JIS Z 2801. ABS interior ladders with the capacity of 9 or 2x9 or 3x9 levels enabling to fill also perforated basket ISO 600x400 as well as with the interior elements of the Eolis range.

Ref. **1007.11 W01** exactly equipped as per picture with:

- 2 small + 2 medium + 1 big drawer
- 4 polypropylene with threadguards Ø 125 mm wheels, 2 with brakes, stainless steel balls bearing.
- Antibacterial top work surface
- 69 W x 49 D x 96,5 H cm



Ref. **1007.12 W01** exactly equipped as per picture with:


- 4 small + 4 medium + 2 big drawers
- 6 polypropylene with threadguards Ø 125 mm wheels, 3 with brakes, stainless steel balls bearing.
- Antibacterial work-top surface
- 132,5 W x 49 D x 96,5 H cm



Ref. **1007.13 W01** exactly equipped as per picture:

- 6 small + 6 medium + 3 big drawers
- 8 polypropylene with threadguards Ø 125 mm wheels, 4 with brakes, stainless steel balls bearing.
- Antibacterial work-top surface
- 196,6 W x 49 D x 96,5 H cm

LOCKING



Key lock

Ref. **1700.01**



Code lock

Ref. **1700.02**



Badge lock

Ref. **1700.03**



Plum



Green



Saffron



Coral



Lilac



Blue



Red

WALL CABINET AND DRAWERS BLOCK

eoli 600 x 400

WALL CABINET



Ref. **1000.41 C01** exactly equipped as per picture with:

- 1 shelf
- 1 curtain with code lock
- overall dimensions: 66.5 W x 43.5 D x 40.5 H cm
- Inside dimensions: 60.8 W x 40 D cm



Other configuration

DRAWERS BLOC



Ref. **1007.22 D01**

Exactly equipped as per picture
• 69 W x 49 D x 73 H cm


- Epoxy steel welded monobloc body frame. Closed on 5 sides - Opening side: 600 mm.
- 4 polypropylene with threadguards swivel Ø 75 mm wheels, 2 with brakes, stainless steel balls bearing.
- ABS interior ladders.



Ref. **1007.20 D01**

Exactly equipped as per picture
• 69 W x 49 D x 85.5 H cm

MOBILE TABLES



Ref. **1000.52 T01** exactly equipped as per picture

- with 3 drawers
- 155W x 57,5D x 91H cm



Ref. **1000.50**

- 155W x 57,5D x 91H cm
- Other size: 200W x 57,5D x 91H cm



Plum



Green



Saffron



Coral



Lilac



Blue



Red

TREATMENT CART WITH ROLLING SHUTTER

Product advantages:

- The accessories can be moved without using tools.
- No metal sliding rails for best hygiene.

eolis 600 x 400

TREATMENT CART WITH BASKETS



Ref: **1006.11 T01** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **key lock**.
- 2 adj. rails left/right + adj. handle.
- 4 wheels Ø 125 mm, 2 brakes.
- 5 baskets (2 small, 2 medium, 1 big) with separations.
- Sticker TREATMENT
- 77W x 59D x 101H

TREATMENT CART WITH BASKETS



Ref: **1006.11 T02** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **key lock**.
- 2 adj. rails left/right + adj. handle.
- 4 wheels Ø 125 mm, 2 brakes.
- 4 baskets (4 medium) with separations.
- Sticker TREATMENT
- 77W x 59D x 101H

TREATMENT CART WITH DRAWERS AND CODE LOCK



Ref: **1006.11 T03** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **electronic lock**.
- 2 adj. rails left/right + adj. handle.
- 4 wheels Ø 125 mm, 2 brakes.
- 5 drawers (2 small, 2 medium, 1 big).
- Sticker TREATMENT
- 77W x 59D x 101H

TREATMENT / MEDICATION CART WITH PATIENT DRAWERS



Ref: **1006.11 T04** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **key lock**.
- 2 adj. rails left/right + adj. handle.
- 4 wheels Ø 125 mm, 2 brakes.
- 15 drawers (12 big, 3 small) with separations.
- Sticker TREATMENT
- 77W x 59D x 101H

ADD EXTRA ACCESSORIES pages 12 to 15



Plum



Green



Saffron



Coral



Lilac



Blue



Red

MEDICATION CART WITH ROLLING SHUTTER

Product advantages:

- The accessories can be moved without using tools.
- No metal sliding rails for best hygiene.

eolis 600 x 400

MEDICATION CART 18 PATIENTS WITH CODE LOCK



Ref: **1006.11 T05** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **electronic lock**.
- 2 adj. rails left/right.
- 4 wheels Ø 125 mm, 2 brakes.
- 18 patient bins + 3 drawers.
- Sticker MEDICATION
- 77W x 59D x 101H

MEDICATION CART 21 PATIENTS WITH KEY LOCK



Ref: **1006.11 T06** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **key lock**.
- 2 adj. rails left/right.
- 4 wheels Ø 125 mm, 2 brakes.
- **21 patient bins** + 3 drawers.
- Sticker MEDICATION
- 77W x 59D x 101H

MEDICATION CART 40 PATIENTS WITH CASSETTE



Ref : **1006.11 T07** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **key lock**.
- 2 adj. rails left/right.
- 4 wheels Ø 125 mm, 2 brakes.
- **40 patient bins** + 1 cassette + 1 basket.
- Sticker TREATMENT
- 77W x 59D x 101H

MEDICATION CART 90 PATIENTS PILLS DISPENSER



Ref : **1006.11 T08** EOLIS cart 600x400 exactly equipped as per picture:

- ABS, top tray with rims.
- Extract. rolling shutter enabling disinfection + **key lock**.
- 2 adj. rails left/right.
- 4 wheels Ø 125 mm, 2 brakes.
- **90 patient** pills dispensers.
- Sticker TREATMENT
- 77W x 59D x 101H

ADD EXTRA ACCESSORIES pages 12 to 15




CREATE YOUR TROLLEY EOLIS 600

1

THE BASE

Without shutter



10 levels



Ref. **1006.11**

10 levels

- Lateral and rear colour ABS panel.
- ABS upper beige tray with rims.
- Rolling shutter with key, code or badge lock in option (closing from top to bottom).
- Extractable rolling shutter enabling the cleaning and the disinfection.
- Integrated height-adjustable push handle included.

2

LOCKING SHUTTER



Rolling shutter only (no lock) ref. **1600.10**

1 level



Rolling shutter + key lock ref. **1600.11**

1 level

3

DRAWERS

1 level - ref. **1600.30**

600 x 400 x 63 mm ht



2 levels - ref. **1600.31**

600 x 400 x 138 mm ht



3 levels - ref. **1600.32**

600 x 400 x 216 mm ht



SEPARATION KITS

Kit for small drawer 1 level

2 cases **1600.70** 8 cases **1600.74**

3 cases **1600.71** 12 cases **1600.75**

4 cases **1600.72** 18 cases **1600.76**

6 cases **1600.73** 24 cases **1600.77**

Kit for medium drawer 2 levels

2 cases **1600.80** 8 cases **1600.84**

3 cases **1600.81** 12 cases **1600.85**

4 cases **1600.82** 18 cases **1600.86**

6 cases **1600.83** 24 cases **1600.87**

Kit for big drawer 3 levels

2 cases **1600.90** 8 cases **1600.94**

3 cases **1600.91** 12 cases **1600.95**

4 cases **1600.92** 18 cases **1600.96**

6 cases **1600.93** 24 cases **1600.97**

SHELVES

Telescopic shelf - ref. **1101.99**

• 1 level



Resin shelf - ref. **1100.58**

• 1 level



Polypropylene shelf - ref. **1100.26**

• With rim - 1 level



BASKETS


1 level ref. **1102.30** - 50 mm ht Separations



2 levels - ref. **1102.31** - 100 mm ht



3 levels - ref. **1102.32** - 200 mm ht



End stop **1102.39** Label holders **1102.69**



Drawers, baskets, patient bins and medication trays can be mixed - maximum 10 levels capacity

4

COLOURS



Plum



Green



Saffron



Coral



Lilac



Blue



Red


5

ADD EXTRA ACCESSORIES (see pages 12 to 15)

x 400 WITH ROLLING SHUTTER


- 2 height-adjustable normalized rails (1 on each side) included.
- 4 swivel Ø 125 mm polypropylene with threadguards wheels, 2 with brakes, stainless steel ball bearing. All around ABS bumpers.
- Weight of an empty cart with rolling shutter: 29.5 kg
- Rolling shutter is taking 1 level, the rests for components.

Ref. **1006.11**
10 levels



THE BASE

Without shutter



1

LOCKING SHUTTER

2



Rolling shutter + Code lock ref. **1600.12**
1 level



Rolling shutter + Badge lock ref. **1600.13**
1 level

PATIENT BINS

MEDICATION TRAYS

3

1 level - ref. **1100.26 OP3**

- 6 patients bins + 6 separations
- 90 x 385 x 60 mm ht each



2 levels - ref. **1100.26 OP6**

- 6 patients bins + 6 separations
- 90 x 385 x 136 mm ht each



1 level - ref. **1100.26 OP4**

- 4 patients bins + 4 separations
- 130 x 385 x 60 mm ht each



2 levels - ref. **1100.26 OP7**

- 4 patients bins + 4 separations
- 130 x 385 x 136 mm ht each



1 level - ref. **1100.26 OP5**

- 3 patients bins + 3 separations
- 178 x 385 x 60 mm ht each



2 levels - ref. **1100.26 OP8**

- 3 patients bins + 3 separations
- 178 x 385 x 136 mm ht each



1 level - ref. **1103.41**

- 7 pills dispenser
- with 3 separations
- 78 x 380 x 55 mm ht



1 level - ref. **1103.38**

- 8 pills dispenser
- with 3 separations
- 65 x 380 x 65 mm ht



1 level - ref. **1103.37**

- 8 pills dispenser
- with 3 separations
- 65 x 380 x 40 mm ht



1 level - ref. **1103.39**

- 7 pills dispenser
- with 4 separations
- 28 x 280 x 40 mm ht



CASSETTE - LABELS

2 levels - ref. **1600.36**

3 levels - ref. **1600.37**



Handle - ref. **1600.42**

Labels - ref. **1103.36**

- Colored sheets of 28 precut labels



Label Holder - ref. **1103.35**

- Pack of 24 label holders



Drawers, baskets, patient bins and medication trays can be mixed - maximum 10 levels capacity

COLOURS

4



Plum



Green



Saffron



Coral



Lilac



Blue



Red

ADD EXTRA ACCESSORIES (see pages 12 to 15)

5

EXTRA ACCESSORIES FOR **eoli** 600 x 400



Sliding lateral shelf - Ref: **1600.43**

• only for trolley with drawer ref 1006.01



Foldaway lateral shelf - Ref: **1500.01**

• For Eolis ref 1006.01 and 1006.11



Telescopic front shelf
Ref: **1600.27**

• For Eolis ref 1006.01 and cabinets



Telescopic front drawer
lockable - Ref: **1600.28**

• For Eolis ref 1006.01 and cabinets



Telescopic front drawer
lockable - Ref: **1600.29**

• For Eolis ref 1006.11 and cabinets



Telescopic epoxy steel shelf
600 x 400 - Ref: **1101.99**

• For Eolis ref 1006.11 and cabinets



Guard rail - Ref: **1600.40**



Guard rail with bottle
holder - Ref: **1600.41**



3 cases bottle holder - Ref: **392.28**



Additional adjustable rail
Ref: **1500.03**

• 330 L x 10 D x 35 H



2 curved hooks adjustable IV pole
Ref: **1500.23**



Universal sharps container holder
Ref: **1101.97**

• without sharp container

EXTRA ACCESSORIES FOR 600 x 400




Gloves box holder in epoxy steel for rail
Ref: **1100.74**



Epoxy steel lateral basket

- 300 L x 100 D x 100 mm ht
Ref: **1101.90**
- 300 L x 160 D x 100 mm ht
Ref: **1101.91**



Side basket in epoxy steel

- 400 L x 200 D x 150 mm ht
Ref: **494.23**
- 400 L x 100 D x 100 mm ht
Ref: **494.21**



Lateral box in epoxy steel for filing pockets - Ref: **1101.95**

- 328 L x 165 D x 245 mm ht



Lateral letter holder A4 size

- 1 case - Ref: **1500.29**
- 2 cases - Ref: **1500.28**
- 3 cases - Ref: **1500.27**



Sheet holder for rail - Ref: **1101.96**



Single waste bin with removable lid and support

- 10 liters - Ref: **1501.24**
- 15 liters - Ref: **1501.25**



Double waste bins with removable lids and support

- 2 x 10 litres - Ref: **1501.26**
- 2 x 15 litres - Ref: **1501.28**



Refuse bin with lid operated by knee
Ref: **1100.80**



Steel refuse bag holder with clamps for standard rail

- 10 to 20 liters - Ref: **1501.20**
- 20 to 30 liters - Ref: **1501.21**



Stickers

- Ref: **1600.44** - Standard sticker: MEDICATION or TREATMENT or ANAESTHESIA
- Ref: **1600.45** - Customized stickers (max. 10 characters)



Additional cost for 4 antistatic wheels with 4 brakes - Ref: **1600.46**

Additional cost for directional wheel with brakes - Ref: **1100.77**

EXTRA ACCESSORIES FOR **eoli** 600 x 400



Boxes holder - Ref: **1600.24**

Bridge - Ref: **1500.11**

Rail for bridge - Ref: **1600.22**

Liquid soap holder - Ref: **1500.21**



Shelf for bridge
Ref: **1600.23**

Line of plastic boxes

- line of 6 boxes
Ref: **1500.16**


- line of 5 boxes
Ref: **1500.17**

- line of 4 boxes
Ref: **1500.18**



Boxes

- 1 litre - 103 L x 165 D x 83 mm HT - Ref: **890.65**
- 3,5 litres - 145 L x 240 D x 125 mm HT - Ref: **890.66**
- 10 litres - 207 L x 345 D x 165 mm HT - Ref: **890.67**



Catheter holder for rail

- 1 storage - Ref: **1101.94**
- 2 storages - Ref: **1101.92**
- 3 storages - Ref: **1101.93**




Swivel lateral defibrillator
shelf - Ref: **1500.75**



Cardiac board with handles - Ref: **1600.25**

- 600 x 430 mm **CE**



Gaz bottle holder Ø 150 mm max.
Ref: **1100.42**



Drawers for toxic drugs with lid and lock

- 1 level (63 mm H) - Ref: **1600.33**
- 2 levels (138 mm H) - Ref: **1600.34**
- 3 levels (216 mm H) - Ref: **1600.35**



Files holder on telescopic rail
Ref: **1102.14**

- Capacity: 11 cases 40 mm
- 4 levels



Files holder with compartments and removable partitions

- For medical files - 18 cases (3 L x 32 D x 36 H) - 5 levels - Ref: **1102.01**
- For X-ray files - 18 cases (3 L x 40 D x 48 H) - 7 levels - Ref: **1102.04**

EXTRA ACCESSORIES FOR **eoli** 600 x 400



Tube Ø 30 mm for computer holder
Ref: **1500.80**

Vesa holder on short pivot - Ref: **1103.55**



Stainless steel shelf for laptop
Ref: **1500.82**

- Height-adjustable and rotating



Fixed arm with universal securized tray for laptop
Ref: **1500.81**

- Height-adjustable and rotating



Arms with universal securized tray for laptop

- Fixed height - Ref: **1103.53**
- Variable height - Ref: **1103.54**



Arms with universal screen holder

- Fixed height - Ref: **1103.56**
- Variable height - Ref: **1103.57**



Compact central unit holder - Ref: **1600.26**



Self-powered cart computer Venus

Dual battery provides 125W power for standing 6 hours of operation run non-stop 7/24 by swapping batteries. Fan-less and no vent hole design reduce the risk of cross-infectial.

- Touch-screen 19" - Ref: **1104.30**
- Touch-screen 22" - Ref: **1104.31**



2 batteries kit
Li-ion 10.8V / 5800mAh
Ref: **1104.34**



Slot battery bank (2/6)

- For 2 batteries 24V / 150W - Ref: **1104.35**
- For 6 batteries 24V / 250W - Ref: **1104.36**



Lithium universal battery (Life P04) 220 V AC - 250W / h
Ref: **1104.52**

Secured support with key lock for battery ref 1104.52
Ref: **1104.51**



Lithium universal battery (Life P04) 220 V AC - 250W / h
Ref: **1504.01**

- Battery integrated in the cart.
- Capacity: 250 W / h for Eolis cart.

CABINET WITH DRAWERS

Product advantages:

- The accessories can be moved without using tools.
- No metal sliding rails for best hygiene.

eoli 600 x 400



Ref: **1007.01C01** Cabinet exactly equipped as per picture:

- Epoxy coated body, 7 colors, ABS interior ladders,
- 12 drawers (3 big, 4 medium, 5 small).
- 1 push bar included.
- Key lock.
- 4 Wheels Ø 125 mm, 2 brakes.
- 68 W x 49 D x 195 HT cm.



Ref: **1007.02C01** - Cabinet exactly equipped as per picture:

- Epoxy coated body, 7 colors, ABS interior ladders,
- 10 drawers (2 big, 4 medium, 4 small).
- 1 push bar included.
- Key lock.
- 4 Wheels Ø 125 mm, 2 brakes.
- 68 W x 49 D x 165 HT cm.

EXTRA ACCESSORIES



Foldable shelf for cabinet
Ref. **1700.24**



Additional push handle
Ref. **1700.21**



Medical rail for cabinet
Ref. **1700.22**



Wheels buffer
Ref. **1700.20**



Plum



Green



Saffron



Coral



Lilac



Blue




Red

CREATE YOUR CABINET WITH DRAWERS

ecoli 600 x 400

1 THE BASE




Ref. **1007.03**

64 x 46 x 134 HT cm
14 levels




Ref. **1007.08**

128 x 46 x 134 HT cm
2 x 14 levels




Ref. **1007.02**

64 x 46 x 164 HT cm
18 levels



Ref. **1007.07**

128 x 46 x 164 HT cm
2 x 18 levels



Ref. **1007.01**

64 x 46 x 195 HT cm
22 levels

Ref. **1007.06**

128 x 46 x 195 HT cm
2 x 22 levels

2 LOCKING



Key lock

Ref. **1700.01**



Code lock


Ref. **1700.02**



Badge lock

Ref. **1700.03**

3 DRAWERS



Drawer 1 level - ref. **1600.30** (without kit)
600 x 400 x 63 mm ht



Drawer 2 levels - ref. **1600.31** (without kit)
600 x 400 x 138 mm ht



Drawer 3 levels - ref. **1600.32** (without kit)
600 x 400 x 216 mm ht

4 SEPARATION KITS

SEPARATION KIT for small drawer 1 level

Kit 2 cases	1600.70	Kit 8 cases	1600.74
Kit 3 cases	1600.71	Kit 12 cases	1600.75
Kit 4 cases	1600.72	Kit 18 cases	1600.76
Kit 6 cases	1600.73	Kit 24 cases	1600.77

SEPARATION KIT for medium drawer 2 levels

Kit 2 cases	1600.80	Kit 8 cases	1600.84
Kit 3 cases	1600.81	Kit 12 cases	1600.85
Kit 4 cases	1600.82	Kit 18 cases	1600.86
Kit 6 cases	1600.83	Kit 24 cases	1600.87

SEPARATION KIT for big drawer 3 levels

Kit 2 cases	1600.90	Kit 8 cases	1600.94
Kit 3 cases	1600.91	Kit 12 cases	1600.95
Kit 4 cases	1600.92	Kit 18 cases	1600.96
Kit 6 cases	1600.93	Kit 24 cases	1600.97

5 COLOURS



Colour choice



Plum



Green



Saffron



Coral



Lilac



Blue



Red

CABINETS WITH ROLLING SHUTTER

Product advantages:

- Easy to clean (infection control).
- All drawers and baskets can be removed and changed without any tool.

eolis 600 x 400

Cabinets format 600 x 400 compatible with all the Eolis range equipments (drawers, baskets, shelves...).

Perfectly smooth surfaces without nook or unevenness, resistant to all disinfectants currently used in hospital.



Ref: **1007.51C01** Cabinet exactly equipped as per picture:

- Epoxy coated body, 7 colors, ABS interior ladders,
- 9 ABS baskets (3 big, 4 medium, 2 small).
- 1 push bar included.
- Extractable rolling shutter enabling perfect disinfection.
- Key lock.
- 4 Wheels Ø 125 mm, 2 brakes.
- Sticker MEDICATION
- 69 W x 57 D x 192 HT cm.

Ref: **1007.56C01** Cabinet exactly equipped as per picture:

- Epoxy coated body, 7 colors, ABS interior ladders,
- 20 drawers (6 big, 10 medium, 2 small).
- 1 push bar included.
- Extractable rolling shutter enabling perfect disinfection.
- Key lock.
- 4 Wheels Ø 125 mm, 2 brakes.
- Sticker MEDICATION
- 134 W x 57 D x 192 HT cm (other heights available: 164 and 134 cm).

For basket separators, go to page 10.



Plum



Green



Saffron



Coral



Lilac



Blue



Red

CABINETS WITH ROLLING SHUTTER

Product advantages:

- Easy to clean (infection control).
- Modularity.
- Strong.

eoli 600 x 400

We have developed a large range of cabinets with different configurations: fully opened or closed, with baskets or drawers. We also offer different pharmacy cabinets.



Ref: **1007.51C02** Medication cabinet exactly equipped as per picture:

- Epoxy coated body, 7 colors, ABS interiors ladders, key lock.
- 38 patient bins + 6 drawers.
- 1 push bar included.
- Extractable rolling shutter enabling perfect disinfection.
- 4 wheels Ø 125 mm, 2 brakes.
- Sticker MEDICATION
- 69 W x 57 D x 192 HT cm (other heights available: 162 and 131 cm).



Ref: **1007.91C01** Medication cabinet exactly equipped as per picture:

- Epoxy coated body, 7 colors, ABS interior ladders, key lock + 2 drawers.
- 192 pill dispensers on 24 trays (8 per tray: Monday, Tuesday... till Sunday).
- 1 push bar included.
- Extractable rolling shutter enabling perfect disinfection.
- 4 wheels Ø 125 mm, 2 brakes.
- Sticker MEDICATION
- 69 W x 57 D x 192 HT cm (other heights available: 162 and 131 cm).



Toxic drugs cabinet Ref: **1101.84**

- With steel door and lock



Probes holder shelf Ref: **1100.25**

- telescopic 60 hooks



Telescopic epoxy steel shelf Ref: **1101.99**



Plum



Green



Saffron



Coral



Lilac



Blue



Red

CREATE YOUR CABINET EOLIS 600


1

THE BASE


- Including rolling shutter



Ref. **1007.53**
69 x 57 x 131 HT cm
13 levels



Ref. **1007.58**
134 x 57 x 131 HT cm
2 x 13 levels



Ref. **1007.52**
69 x 57 x 162 HT cm
17 levels



Ref. **1007.57**
134 x 57 x 162 HT cm
2 x 17 levels

2

LOCKING



Key lock - ref. **1700.10**



Code lock - ref. **1700.11**



Seal lock - ref. **1700.12**

3

DRAWERS

1 level - ref. **1600.30**
600 x 400 x 63 mm ht



2 levels - ref. **1600.31**
600 x 400 x 138 mm ht



3 levels - ref. **1600.32**
600 x 400 x 216 mm ht



SEPARATION KITS

Kit for small drawer 1 level

2 cases 1600.70	8 cases 1600.74
3 cases 1600.71	12 cases 1600.75
4 cases 1600.72	18 cases 1600.76
6 cases 1600.73	24 cases 1600.77

Kit for medium drawer 2 levels

2 cases 1600.80	8 cases 1600.84
3 cases 1600.81	12 cases 1600.85
4 cases 1600.82	18 cases 1600.86
6 cases 1600.83	24 cases 1600.87

Kit for big drawer 3 levels

2 cases 1600.90	8 cases 1600.94
3 cases 1600.91	12 cases 1600.95
4 cases 1600.92	18 cases 1600.96
6 cases 1600.93	24 cases 1600.97

SHELVES

Telescopic shelf - ref. **1101.99**

- 1 level



Resin shelf - ref. **1100.58**

- 1 level



Polypropylene shelf - ref. **1100.26**

- With rim - 1 level



BASKETS

1 level ref. **1102.30** - 50 mm ht Separations



Length **1102.33**
Width **1102.34**

2 levels - ref. **1102.31** - 100 mm ht



Length **1102.35**
Width **1102.36**

3 levels - ref. **1102.32** - 200 mm ht



Length **1102.37**
Width **1102.38**

End stop **1102.39** Label holders **1102.69**



Drawers, baskets, patient bins and medication trays can be mixed

4

COLOURS



Plum



Green



Saffron



Coral



Lilac



Blue




Red


5

ADD EXTRA ACCESSORIES (see pages 12 to 15)

x 400 WITH ROLLING SHUTTERS



Ref. **1007.51**
69 x 57 x 192 HT cm
21 levels



Ref. **1007.56**
134 x 57 x 192 HT cm
2 x 21 levels

High capacity thin scale cabinet **only** for pills dispenser 1103.37 and 1103.39 (includ. rolling shutters).

• Including rolling shutter



Ref. **1007.93**
18 levels

1007.98
2 x 18 levels

1007.92
24 levels

1007.97
2 x 24 levels

1007.91
30 levels

1007.96
2 x 30 levels



Key lock - ref. **1700.10**



Code lock - ref. **1700.11**



Seal lock - ref. **1700.12**

PATIENT BINS

1 level - ref. **1100.26 OP3**

- 6 patients bins + 6 separations
- 90 x 385 x 60 mm ht each



2 levels - ref. **1100.26 OP6**

- 6 patients bins + 6 separations
- 90 x 385 x 136 mm ht each



1 level - ref. **1100.26 OP4**

- 4 patients bins + 4 separations
- 130 x 385 x 60 mm ht each



2 levels - ref. **1100.26 OP7**

- 4 patients bins + 4 separations
- 130 x 385 x 136 mm ht each



1 level - ref. **1100.26 OP5**

- 3 patients bins + 3 separations
- 178 x 385 x 60 mm ht each



2 levels - ref. **1100.26 OP8**

- 3 patients bins + 3 separations
- 178 x 385 x 136 mm ht each



MEDICATION TRAYS

1 level - ref. **1103.41**

- 7 pills dispenser
- with 3 separations
- 78 x 380 x 55 mm ht



1 level - ref. **1103.38**

- 8 pills dispenser
- with 3 separations
- 65 x 380 x 65 mm ht



1 level - ref. **1103.37**

- 8 pills dispenser
- with 3 separations
- 65 x 380 x 40 mm ht



1 level - ref. **1103.39**

- 7 pills dispenser
- with 4 separations
- 28 x 280 x 40 mm ht



CASSETTE - LABELS

2 levels - ref. **1600.36**

3 levels - ref. **1600.37**



Handle - ref. **1600.42**

Labels - ref. **1103.36**

- Colored sheets of 28 precut labels



Label Holder - ref. **1103.35**

- Pack of 24 label holders



Drawers, baskets, patient bins and medication trays can be mixed

COLOURS



Plum



Green



Saffron



Coral



Lilac



Blue



Red

ADD EXTRA ACCESSORIES (see pages 12 to 15)

TRANSFER CABINETS



Ref: **1002.26T01** equipped exactly as per picture.

- 12 ABS baskets ISO standard 600x400 (5 small, 4 medium, 3 large) + labels. Rails with end stops. 2 pushbars. Wheels Ø 125 mm.
- 74 W x 50 D x 194 H.



Ref: **1002.29T01** equipped exactly as per picture.

- 16 ABS baskets ISO standard 600x400 (2 small, 6 medium, 8 large) + labels. 2 push bars. Rails with end stops. 4 wheels Ø 125 mm.
- 97 W x 70 D x 194 H.



Ref: **1002.04D01** equipped exactly as per picture.

- 5 ABS baskets ISO standard 600x400 (3 small, 1 medium, 1 big) + dividers for top tray
- 1 push bar. Adjustable rails with end stops. 4 Wheels Ø 125 mm.
- 70 W x 51 D x 97 H.



Ref: **1002.62T01** equipped exactly as per picture.

- Strong aluminium structure. 2 push bars. Doors with key lock opening 270°.
- 20 ABS baskets ISO standard 600x400 (10 small, 6 medium, 4 big). Rails with end stops. 4 wheels Ø 200 mm, 2 with brakes.
- 142 W x 76 D x 179 H

COMPUTER CART

Product advantages:

We validate with you your computer cart project (software compatibility, battery, wireless...).
Tel: 0033 2 43 82 95 69 for assistance.

Our computer work stations are height adjustable, easily mobile and ready to receive your laptop or screen computer.
A battery as option optimize the autonomy.



Cart for laptop (without computer).

- Fixed height computer cart for laptop (105 cm) **Ref 340.01**
- Height adjustable computer cart for laptop (77 to 106 cm HT) **Ref 340.02**
- Secure the computer without tools
- Wheels Ø 100 mm, 2 with brakes.
- Top tray 55 x 55 cm.




Computer cart with vesa holder for screen (without computer and battery).

- Fixed height cart for screen (105 cm) **Ref 340.03**
- Height adjustable cart for screen (77 to 106 cm HT) **Ref 340.04**
- Wheels Ø 100 mm, 2 with brakes.
- Top tray 55 x 55 cm.

EXTRA ACCESSORIES

Waterproof keyboard with touch-pad

1104.16



Drawer + key lock

340.50



File holder

1101.96



Basket 40 x 11 x 7,5 HT cm

494.21

