

Optically clear window for scanning and viewing contents

Optically clear window

Caregivers can now see with clarity the text on the fluid bag and the fluid level in the bag. Can also scan the fluid contents bar code.

Color coded pressure gauge - with 360° visibility

InfusaScan™ incorporates a color coded pressure gauge that can be easily seen from a distance.

Patient dedicated - to help reduce hospital infections

InfusaScan has a place for the patient's name, date, and fluid contents. Single patient use may reduce the potential for cross contamination.

Large inflation bulb – for rapid inflation

InfusaScan uses a large capacity, rapid recovery inflation bulb, which shortens inflation time and reduces hand fatigue.

Latex and neoprene free

Patients and clinicians are protected from latex/neoprene allergies.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.


GE, the GE Monogram, Vital Signs and InfusaScan are trademarks of General Electric Company.

Vital Signs, Inc.
A GE Healthcare Company
20 Campus Road
Totowa, NJ 07512
USA


www.vital-signs.com

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
Tel. +358 10 394 11
Fax +358 9 146 3310

www.gehealthcare.com



Infusascan's clear window



Netting backed infusor

InfusaScan™ Single Patient Use Pressure Infusor Ordering Information

Product number	Description	Volume (mL)	Qty per Case
IN800020	InfusaScan™ Pressure Infusor	500	20
IN900020	InfusaScan Pressure Infusor	1,000	20



GE imagination at work

M1177176 04/09 © 2009 All rights reserved • U.S. Patent Pending Vital Signs, Inc., a General Electric company, doing business as GE Healthcare. TD152 Rev 01 030209