

Infusable®

Pressure Infusor

- Patient-dedicated to minimize the spread of nosocomial infections via touch contamination.
- Large inflation bulb has rapid recovery for quick and easy inflation.
- Easy-to-use stopcock and bulb-to-tube fitted seal prevent leaks.
- High visibility patented gauge accurately indicates bladder pressure.
- Nylon netting makes loading easy and permits a clear view of fluid level



Product Description	Single patient use pressure infusor		
Ordering Information	Product Code	Description	Qty/Case
	IN800012	500 mL Infusable	12
	IN800048	500 mL Infusable	48
	IN900012	1000 mL Infusable	12
	IN900048	1000 mL Infusable	48
	IN950012	3000 mL Infusable	12



About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our “healthymagination” vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality and efficiency around the world.

Vital Signs, Inc.
A GE Healthcare Company
20 Campus Road
Totowa, NJ 07512
USA

GE Healthcare
P.O. Box 900, FIN-00031 GE, Finland
Tel. +358 10 394 11
Fax +358 9 146 3310

www.gehealthcare.com

www.vital-signs.com

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, the GE Monogram, Vital Signs and Infusable are trademarks of General Electric Company.

M1198767 01/10 © 2010 All rights reserved
TD058 Rev.02 051704 © 2004 Vital Signs, Inc.
Patent No. 4,735,613 Vital Signs, Inc. a General Electric company, doing business as GE Healthcare.



GE imagination at work