


Vital Signs
a GE Healthcare Company



GreenLight II* Laryngoscope with
Single-Use Steelite Blades



Vital Signs Laryngoscopes

Designed to Reduce Infection

- Strength of steel, safety of single use
- Single-use blades help reduce the risk of infection and cross-contamination
- Compatible with Safe Sac*, a disposable handle sheath that helps keep handle clean



Safe Sac Handle Sheath



Traditional Handle




Stubby Handle

Designed for Performance

- Ergonomic, over-molded handles are lightweight and designed to fit the clinician's hand
- Handle has a durable, slip-resistant grip and a wide, sturdy base
- Bright LED light with a covered light rod prevents clinician back-flash glare

Designed to make things Easy

- Single-use blades reduce time and effort associated with reprocessing of reusable blades
- Standard AA batteries (Alkaline or Lithium) used in traditional handle for easy replacement (Single CR-123 Lithium battery used in stubby handle)
- Never worry about replacing bulbs again with our laryngoscope design



Steelite Blades (Plastic Blades also available)

Designed to Last

- Efficient technology extends battery life beyond that of a standard incandescent bulb light source
- Five-year warranty on all handles

Designed for You

ABOUT GE HEALTHCARE:

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our “healthymagination” vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world. Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

Vital Signs, Inc.
a GE Healthcare Company
20 Campus Road
Totowa, NJ 07512
USA

© 2012 General Electric Company – All rights reserved. General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram, and Vital Signs are trademarks of General Electric Company.

*GreenLight II and Safe Sac are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.

L2008_GreenLight Brochure_EN_v4_July 2012

DOC1075612



imagination at work