

Laryngoscope Systems

Greenlight II

ISO Green Standard
Laryngoscope Handle and Blades

Vital View II

Laryngoscope Handle
and Blades



The Greenlight II™ and Vital View™ II provide a wide range of features and benefits for laryngoscopy which may enhance clinical application, patient care, and materials management. Both handles are backed by a five year warranty and are compatible with the Safe Sac™ laryngoscope handle cover, which helps to prevent soiling and cross-contamination.

Bright: The LED (Light Emitting Diode) in the Greenlight II and Vital View II is one of the brightest light sources on the market. The proximal end of the light rod on each laryngoscope blade is covered to prevent reflected light, which may interfere with visualization.

Long-Lasting: The LED in both systems lasts approximately 50,000 hours. The efficiency of the Greenlight II and Vital View II technology extends battery life up to seven times longer than standard incandescent bulb light sources. Also, both systems provide constant light output with battery decline, whereas traditional laryngoscope systems may decrease in illumination over time.

Ergonomic: The over-molded handles of both laryngoscope systems have a durable, slip resistant grip.

Quality of Care: Helps to mitigate the risk of cross-contamination and infection by using Vital Signs single use stainless steel or engineered plastic laryngoscope blades. They provide consistent and optimal illumination. Reusable blades may decrease in light output due to degradation from multiple sterilization cycles.

Productivity: Single patient use blades may enhance productivity and reduce costs associated with the several steps of decontamination, sterilization and transportation to reprocess reusable blades. Also, single use blades eliminate the need to replace lost and/or worn reusable blades.



Greenlight II™ Laryngoscope Handle and Blades (ISO Green Standard)



VitalView™ II Laryngoscope Handle and Blades

	Descriptions	Case Qty.	Product Codes
Greenlight II™ Laryngoscope System* (ISO Standard)†	Standard Handle	1	4558GSP
	Stubby Handle	1	4559GSP
	Plastic Macintosh Blades sizes 2 / 3 / 4	20	4002GSP / 4003GSP / 4004GSP
	Plastic Miller Blades sizes 0 / 1 / 2 / 3	20	4010GSP / 4011GSP / 4012GSP / 4013GSP
	Stainless Steel Macintosh Blades sizes 2 / 3 / 4	20	4602 / 4603 / 4604
	Stainless Steel Miller Blades sizes 0 / 1 / 2 / 3	20	4610 / 4611 / 4612 / 4613
Vital View™ II Laryngoscope System‡	Standard Handle	1	4556
	Plastic European Macintosh Blades sizes 2 / 3 / 4	20	4002 / 4003 / 4004
	Plastic Miller Blades sizes 0 / 1 / 2 / 3	20	4510 / 4511 / 4512 / 4513
	Plastic High Visibility Macintosh Blades sizes 2 / 3 / 4	20	4502 / 4503 / 4504
Safe Sac™ Handle Cover	Laryngoscope Handle Cover	400	4552

Vital Signs, Inc.
A GE Healthcare Company
20 Campus Road
Totowa, NJ 07512
USA
www.vital-signs.com

GE Healthcare
P.O. Box 900, FIN-00031 GE,
Finland
Tel. +358 10 394 11
Fax +358 9 146 3310
www.gehealthcare.com

* The Greenlight II handle can be steam autoclaved (please refer to the handle's Instructions for Use).
† Technical Committee ISO/TC 121, A. a. (2003). ISO STANDARD 7376:2003(E) Anaesthetic and respiratory equipment - Laryngoscopes for tracheal intubation. Switzerland: International Organization for Standardization.
‡ Vital View II is a non-standard handle and blade configuration and is only available with plastic laryngoscope blades.

© 2009 General Electric Company - All rights reserved. General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.
GE, GE Monogram, Vital Signs, Vital View, Safe Sac and Greenlight II are trademarks of General Electric Company.
Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.
EMEA M1193253, based on the global code L2008
Greenlight Brochure A4 EN Rev. Rel. 10/2009



GE imagination at work