

Vital Signs  
A GE Healthcare Company

# Limb-O™

When less  
is more



The innovative streamlined design  
of the Limb-O makes it the ideal  
breathing circuit for ventilation.



- Weight** Limb-O weighs 25% less than traditional two-limb circuits which reduces torque on the patient's airway.
- Inventory** Limb-O consolidates your inventory! The size of the lumen accommodates both adult and pediatric patients.
- Compliant** Limb-O's low compliance tubing allows for delivery of a more accurate tidal volume.
- Heat Loss** Limb-O's common wall/septum separating inspiratory and expiratory lumens permits heat transfer from warm exhaled gasses. The innovative streamlined design of the Limb-O makes it the ideal breathing circuit for ventilation.

Limb-O offers versatility with this single limb circuit when used as a transport circuit and also in the OR.

Product Description		<b>Limb-O™ Adult and Pediatric Volume Ventilator Circuit</b>											
Ordering Information (20/case)		Product Code				Description (Code assignment depends on configuration; basic configurations are as follows)							
		AFNXXXX X				1.8 m length							
Accessories		5510				Limb-O Transport System - 7' Star Lumen Oxygen Tubing set and APL Valve Set							
Size		1.8 m											
Weight with Connectors		170 g											
Compliance of Circuit @ 100 cm H <sub>2</sub> O <sup>1</sup>		.93 ml/cmH <sub>2</sub> O											
Insp. & Exp. Resistance at Various Flow Rates <sup>1</sup> (cmH <sub>2</sub> O)	Length	10 L/min		20 L/min		30 L/min		40 L/min		50 L/min		60 L/min	
		I	E	I	E	I	E	I	E	I	E	I	E
	1.8 m	.14	.16	.46	.57	1.06	1.24	1.80	2.12	2.74	3.28	3.69	4.39
Circuit Materials		Corrugated Tubing and Septal Wall - Polyethylene/EVA 22 mm Connectors - Medical Grade Polypropylene Cuff Band - Polyisoprene Expandable Extension Hose - Polypropylene											

Vital Signs, Inc.  
A GE Healthcare Company  
20 Campus Road  
Totowa, NJ 07512  
USA

[www.vital-signs.com](http://www.vital-signs.com)

GE Healthcare  
P.O. Box 900, FIN-00031 GE, Finland  
Tel. +358 10 394 11  
Fax +358 9 146 3310

[www.gehealthcare.com](http://www.gehealthcare.com)

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, the GE Monogram, Vital Signs and Limbo-O are trademarks of General Electric Company.

M1182298 06/09 © 2009 All rights reserved • U.S. Patent Pending Vital Signs, Inc., a General Electric company, doing business as GE Healthcare. TD075 Rev.01 021203



GE imagination at work