


Bag-extractor Bio15



Bag-extractor Bio 15

Bag-extractor Bio15 is intended to make easy blood component separation in any condition.

Bag-extractor Bio15 is suitable both for automatic routine procedures and repeated manual operations, giving a wide and safe component separation and the operator can easily pull the bag set apart.

Bag-extractor Bio15 is the new generation of manual plasma extractor providing the manual-controlled separation of erythrocytes (RCC), buffy-coat (BC), platelet poor plasma (PPP) or platelet rich plasma (PRP) and platelet concentrates.

Technical Data

Max volume:	800 mL.
Dimensions:	L 160 x W 270 x H 305
Weight:	4 Kg
Classification:	Class I