

Vital Signs
a GE Healthcare Company

RediTube*

Tracheal Tube with Pre-Loaded
and Sized Stylet



RediTube: Tracheal Tube with Pre-Loaded and Sized Stylet

RediTube helps ensure you are prepared when seconds count. Integrating the tracheal tube and stylet in one product means less set up time and more convenience. Furthermore, RediTube can reduce the expense and effort of purchasing stylets and tracheal tubes separately.

RediTube Tracheal Tube

- High volume, low pressure style cuff
- Murphy eye
- Wide range of sizes: cuffed and uncuffed

RediTube Stylet

- Pre-loaded and sized for proper position in tracheal tube
- Malleable for adjusting shape of tracheal tube
- Helps support tracheal tube placement

Convenience, Confidence and Cost Effectiveness

Stylets and tracheal tubes are also available separately from Vital Signs.

Vital Signs, Inc.
a GE Healthcare Company
20 Campus Road
Totowa, NJ 07512
USA

© 2012 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, the GE Monogram and Vital Signs are trademarks of General Electric Company.

*RediTube is a trademark of General Electric Company

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.

L2001 RediTube_Broch_EN_Global_v1_January 2012



imagination at work

DOC1086735