

LipiGuard® SB Reinfusion Filter for Salvaged Blood



Reduces fat, leucocyte, anaphylatoxin C3a and microaggregate debris found in salvaged blood

- Combines LipiGuard filter media technology with 40 µm lock woven microaggregate screen
- Designed for use with either intra- or post-operatively salvaged blood that may be washed or unwashed, thereby reducing inventory requirements for a number of different filter products
- Low volume (20 mL) housing allows excellent blood product recovery
- High flow rate for rapid infusion
- Easy spike insertion with a secure hold and rapid filter priming
- Effectively replaces conventional microaggregate filters
- Latex free