

Syringe Driver

SD-10/ 99v1



★ **Selectable infusion mode:**

1- Variable infusion speed: 1-10 or 1-99mm/hr

2- Syringe fully infused in 6, 8, 10 or 12hrs

★ **Ultra low power consumption**

★ **Charging battery from any USB port**

★ **User friendly keypad**

★ **Bolus function**

★ **Audible and visual alarms**

 1254



ISO 9001

ISO 13485

FDA QSR 21

CFR 820 cGMP

Ultra Low Power Consumption