

DIMAR

Medical Devices

DimAir Line

Ambulance
Emergency
Department



CPAP COMFORT COULISSE

AMBULANCE, FIRST AID, EMERGENCY

Single patient Helmet for the treatment of acute respiratory diseases in C.P.A.P. in Emergency.

As it is provided with fast “Coulisse” closing system, it can be quickly placed by one operator only.

Being One size, it ensures the pneumatic seal on patients of any size.

With its reduced “Pocket” packaging, it takes up a little space in Ambulance.

Provided with highly efficient “Easy Vent” Venturi generator, adjustable PEEP valve and manometer to check the pressure, it is directly connected to the Oxygen cylinder, to treat patients with Cardiogenic Aedema in C.P.A.P., starting from their house.

Excellent relationship quality/price.

TECHNICAL FEATURES

- One size
- “Coulisse” closing system
- Venturi generator, PEEP valve and manometer already connected for a prompt use
- Soft straps with quick fastening system
- Made of biocompatible plastic material
- Single patient use
- Latex Free


ADVANTAGES

- Quick placement, opening and closing operations
- It can be placed by one operator only
- Quick airways pressurization
- It can be applied on patients of any size
- Minimal amount of space “Pocket” packaging, which can be easily placed everywhere
- Effective and practical
- Excellent relationship quality/price

CODE: 500/8612

- Single pack
- Multiple pack with 10 pieces

Possibility to customize the different models, versions and packaging in response to customers' specific requirements.



DIMAR

Via G. Galilei, 6 41036 Medolla -MO- I
Tel: +39 0535 611336 Fax: +39 0535 611328
Mail: dimar.srl@tiscali.it

